



Specification

1.00468.0100 L(+)-Ascorbic Acid for analysis EMSURE® ACS, Reag. Ph Eur

Specification		
Assay (iodometric)	99.0 - 100.5	%
Identity (IR-spectrum)	conforms	
Appearance	white or almost white, crystalline powder	
Appearance of solution (50 g/l CO ₂ -free water)	clear (≤ 3 NTU) and not so intense in colour than reference solution BY ₇	
pH (50 g/l CO ₂ -free water)	2.1 - 2.6	
Spec. rotation [α] _D ²⁰ (100 g/l, water)	+20.5 - +21.5	°
Chloride (Cl)	≤ 50	ppm
Sulfate (SO ₄)	≤ 20	ppm
Cu (Copper)	≤ 5	ppm
Fe (Iron)	≤ 2	ppm
Heavy metals (ACS)	≤ 10	ppm
Oxalic acid	≤ 0.2	%
Related substances (HPLC) (Impurity C)	≤ 0.15	%
Related substances (HPLC) (Impurity D)	≤ 0.15	%
Related substances (HPLC) (unspecified impurities singly)	≤ 0.10	%
Related substances (HPLC) (sum of impurities (except impurity C and D))	≤ 0.2	%
Sulfated ash (600 °C)	≤ 0.05	%
Loss on Drying (105 °C)	≤ 0.1	%

Dr. Sebastian Lips

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.