



Specification

1.00807.1000 Trichloroacetic acid for analysis EMSURE® ACS, Reag. Ph Eur

Specification		
Assay (alkalimetric)	≥ 99.5	%
Identity (IR-spectrum)	passes test	
Appearance	colorless, deliquescent crystals	
Appearance of solution (200 g/l; water)	clear and colorless	
In water insoluble matter	≤ 0.01	%
Chloride (Cl)	≤ 10	ppm
Nitrate (NO ₃)	≤ 20	ppm
Phosphate (PO ₄)	≤ 5	ppm
Sulphate (SO ₄)	≤ 200	ppm
Heavy metals (as Pb)	≤ 20	ppm
Cu (Copper)	≤ 5	ppm
Fe (Iron)	≤ 10	ppm
Readily carbonisable substance	passes test	
Sulfated ash (600 °C)	≤ 300	ppm

Dr. Ralf Burgert

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.