



# Specification

---

1.01988.1000 1-Butanol for liquid chromatography LiChrosolv®

---

	Specification	
Purity (GC)	≥ 99.8	%
Identity (IR)	conforms	
Evaporation residue	≤ 2.0	mg/l
Water	≤ 0.05	%
Acidity	≤ 0.0002	meq/g
Alkalinity	≤ 0.0002	meq/g
Transmission (at 230 nm)	≥ 75	%
Transmission (at 240 nm)	≥ 85	%
Transmission (from 310 nm)	≥ 99	%

Filtered by 0.2 µm filter

Dr. Ute Volkwein

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.