



Specification

1.02038.0100 Cesium chloride for analysis EMSURE®

Specification Values

Assay (argentometric)	min 99.5	%
Sulphate (SO ₄)	max 0.002	%
Total nitrogen (N)	max 0.001	%
Al (Aluminium)	max 0.0005	%
Cu (Copper)	max 0.0003	%
Fe (Iron)	max 0.0003	%
K (Potassium)	max 0.002	%
Li (Lithium)	max 0.00005	%
Mg (Magnesium)	max 0.0005	%
Na (Sodium)	max 0.002	%
Pb (Lead)	max 0.0001	%
Rb (Rubidium)	max 0.008	%
Zn (Zinc)	max 0.0003	%

Dr. Andreas Lang

responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature