



# Specification

---

1.02378.0500 Calcium chloride anhydrous granules Reag. Ph Eur

---

	Spec. Values	
Assay (calculated on dried substance)	≥ 98.0	%
Free alkali (as Ca(OH) <sub>2</sub> )	≤ 0.2	%
Water absorption capacity (24 hrs., 80 % relative humidity)	≥ 25.0	%
Loss on drying (200 °C)	≤ 5.0	%

Conforms Reag. Ph.Eur.

Claudia Wiegand  
Responsible laboratory manager quality control

---

This document has been produced electronically and is valid without a signature.