



Specification

1.02790.0250 Copper(II) sulfate pentahydrate for analysis EMSURE® ACS,ISO,
Reag. Ph Eur

	Spec. Values	
Assay (iodometric)	99.0 - 100.5	%
Insoluble matter	≤ 0.005	%
Chloride (Cl)	≤ 0.0005	%
Total nitrogen (N)	≤ 0.001	%
Ca (Calcium)	≤ 0.005	%
Fe (Iron)	≤ 0.003	%
K (Potassium)	≤ 0.001	%
Na (Sodium)	≤ 0.005	%
Ni (Nickel)	≤ 0.005	%
Pb (Lead)	≤ 0.005	%
Zn (Zinc)	≤ 0.03	%

Corresponds to ACS,ISO,Reag. Ph Eur

Claudia Wiegand

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.