



Specification

1.03098.0005 2,2'-Bipyridine GR for analysis (reagent for iron(II) and molybdenum) ACS

	Specification	
Assay (perchloric acid titration)	≥ 99.5	%
Identity (IR-spectrum)	passes test	
Melting point	70 - 72	°C
Absorbance of solution (50% w/v; methanol)	conforms	
Sensitivity (to ferrous ion)	passes test	

Ioannis Chartomatsidis

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.