



Specification

1.04391.1000 n-Hexane for liquid chromatography LiChrosolv®

	Specification	
Purity (GC)	≥ 98.0	%
Identity (IR)	conforms	
Evaporation residue	≤ 1.0	mg/l
Water	≤ 0.01	%
Acidity	≤ 0.0002	meq/g
Alkalinity	≤ 0.0002	meq/g
Transmission (at 210 nm)	≥ 50	%
Transmission (at 220 nm)	≥ 85	%
Transmission (from 245 nm)	≥ 98	%

Filtered by 0.2 µm filter

Jeannette David

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.