



Specification

1.05426.1011 Ammonia solution 32% EMPLURA®

	Specification	
Assay (acidimetric, NH ₃)	≥ 30	%
Identity	passes test	
Color	≤ 10	Hazen
Chloride (Cl)	≤ 0.002	%
Sulfate (SO ₄)	≤ 0.003	%
Heavy metals (as Pb)	≤ 0.0001	%
Fe (Iron)	≤ 0.0005	%
Substances reducing potassium permanganate (as O)	≤ 0.0008	%
Residue on ignition (as sulfates)	≤ 0.003	%

Dr. Reiner Vogt

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.