



# Specification

---

1.06580.0500 di-Sodium hydrogen phosphate dihydrate for analysis EMSURE®

---

	Spec. Values	
Assay (acidimetric)	≥ 99.5	%
Identity	passes test	
pH-value (5 %; water)	9.0 - 9.2	
Chloride (Cl)	≤ 0.001	%
Sulfate (SO <sub>4</sub> )	≤ 0.005	%
Total nitrogen (N)	≤ 0.001	%
Heavy metals (as Pb)	≤ 0.001	%
Cu (Copper)	≤ 0.0003	%
Fe (Iron)	≤ 0.001	%
K (Potassium)	≤ 0.005	%
Pb (Lead)	≤ 0.001	%

Claudia Wiegand

Responsible laboratory manager quality control

---

This document has been produced electronically and is valid without a signature.