



Specification

8.00947.1000 Ethylenediamine for synthesis

Specification

Color according to color reference solution Ph.Eur.	colorless liquid	
Assay (GC, area%)	≥ 99.0	% (a/a)
Water (GC, area%)	≤ 1.0	% (a/a)
Density (d 20 °C/ 4 °C)	0.896 - 0.898	
Identity (IR)	passes test	

Due to its specific melting range the product may be solid, liquid, a solidified melt or a supercooled melt.

Dr. Jörg Bauer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.