

AUSTRALIAN CHEMICAL REAGENTS  
**SAFETY DATA SHEET**

Date Prepared: June 2017  
Version No: 5

---

## 1. IDENTIFICATION OF THE MATERIAL AND SUPPLIER

---

Product Name: AUSDET WBC  
Product Code: 0801  
Other Names: Nil  
Uses: Used in water baths to control algae and slimes. Product is diluted 1 part to 200 for use.

Supplier: Australian Chemical Reagents  
38-50 Bedford Street Gillman SA 5013

Contacts: Telephone: 61 08 84402000  
Fax: 61 08 84402001  
Emergency Phone: 61 08 84402000 Mon-Fri 8:30am - 5:00pm

---

## 2. HAZARDS INFORMATION

---

**Hazard classification:** Classified as non-Hazardous according to the Globally Harmonised System of classification and labelling of Chemicals (GHS) including Work, Health and Safety regulations, Australia.

---

## 3. COMPOSITION / INFORMATION ON INGREDIENTS

---

### Ingredients :

Chemical Entity	CAS No	Proportion
Benzylkonium chloride	[ 68424-85-1 ]	<5%
Water	[7732-18-5]	to 100%

---

## 4. FIRST AID MEASURES

---

Safety showers and eye wash facilities should be provided.

### **Swallowed :**

If conscious wash out mouth with water. Seek medical advice. Show this SDS to medical practitioner.

### **Eye :**

Immediately hold eyelids open and flood with water for at least 15 minutes. Obtain medical aid. Show this SDS to medical practitioner.

### **Skin :**

Remove contaminated clothing. Immediately wash skin thoroughly with water and mild soap. Seek medical advice if irritation persists. Show this SDS to medical practitioner. Launder clothing before reuse.

### **Inhaled :**

Remove from contaminated air. Maintain breathing with artificial respiration if necessary. Seek medical assistance. Show this SDS to a doctor.

---

## 5. FIRE FIGHTING MEASURES

---

### **Suitable Extinguishing Media:**

Water spray. Carbon dioxide, dry chemical powder, or appropriate foam.

### **Hazards From Combustion Products:**

Decomposition products include oxides of carbon.

### **Precautions For Fire Fighters and Special Protective Equipment:**

Fire fighters and others who may be exposed to combustion products during fire should wear full protective clothing including positive pressure self-contained breathing apparatus (SCBA). Wear SCBA with full face-piece, operated in positive pressure mode when fighting fires.

---

## 6. ACCIDENTAL RELEASE MEASURES

---

### Emergency procedures:

Prevent from entering waterways. Restrict access to area. Ventilate area. Remove chemicals that can react with the spilled material.

### Methods and materials for containment and clean up:

Use inert material such as sand or earth to contain spill or leak. Absorb spills with chemical absorber or vermiculite and dispose of in accordance with local regulations.

---

## 7. HANDLING AND STORAGE

---

### Precautions for Safe Handling:

Do not get in eyes, on skin, on clothing. Avoid prolonged or repeated exposure.

### Conditions for Safe Storage:

Store sealed in original container in a cool well ventilated situation away from foods and other chemicals. Do not store in direct sunlight. Observe good hygiene and housekeeping practices.

---

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

---

### National Exposure Standards:

SWA – None known

**Biological Limit Values:** No data available.

### Engineering Controls:

Not required with normal use.

### Personal Protective Equipment (PPE):

The use of rubber gloves complying with AS 2161 and the use of faceshield, chemical goggles or safety glasses with side shield protection complying with AS/NZS 1337 is recommended.

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

---

<b>Appearance :</b>	Clear liquid
<b>Odour:</b>	Nil
<b>pH:</b>	Not known
<b>Boiling Point (°C) :</b>	100
<b>Freezing/melting Point:</b>	Not applicable
<b>Vapour Pressure (mm of Hg @ 25°C) :</b>	Not applicable
<b>Vapour Density:</b>	Not applicable
<b>Specific Gravity :</b>	1
<b>Flash Point (°C) :</b>	Not flammable
<b>Flammability Limits (%) :</b>	Not flammable
<b>Solubility in Water (g/L) :</b>	Soluble

---

## 10. STABILITY AND REACTIVITY

---

### Chemical stability:

Stable.

### Conditions to avoid:

Excessive heat. Sunlight

### Incompatible materials:

Nil

**Hazardous decomposition products:**

Refer to section 5 (Fire Fighting Measures).

**Hazardous reactions:**

Hazardous polymerization will not occur.

---

**11. TOXICOLOGICAL INFORMATION**

---

**Health Effects:**

**Swallowed** : May cause burning taste. May lead to gastric irritation and nausea. For benzylkonium chloride LD50 oral-rat 240mg/kg.

**Eye** : Irritating to eye tissue.

**Skin** : May irritate skin tissue. May defat skin.

**Inhaled** : Vapour may result in dizziness, headaches and nausea.

**Chronic Effects:** Repeated skin contact may result in dermatitis.

---

**12. ECOLOGICAL INFORMATION**

---

**Ecotoxicity:**

No data available.

**Persistence and degradability:**

No data available.

**Mobility:**

No data available.

---

**13. DISPOSAL CONSIDERATIONS**

---

Contact a licensed professional waste disposal service to dispose of this material. Observe all federal, state and local environmental regulations.

---

**14. TRANSPORT INFORMATION**

---

**UN Number:** Not applicable

**UN Proper Shipping Name:** Not applicable

**Class and subsidiary risk(s):** Not applicable

**Packing Group:** Not applicable

**Hazchem Code:** Not applicable

**Special precautions for user** : Nil

---

**15. REGULATORY INFORMATION**

---

**Standard for the Uniform Scheduling of Drugs and Poisons (SUSDP):**

Nil

---

**16. OTHER INFORMATION**

---

**Disclaimer:**

All information given by the Company is offered in good faith and is believed to the best of our knowledge to be accurate. However this information is offered without warranty representation inducement or licence and the Company does not assume legal responsibility for reliance upon the same.

Every person dealing with the materials referred to herein does so at his or her own risk absolutely and must make independent determinations of suitability and completeness of information from all sources to ensure their proper use.