

AUSTRALIAN CHEMICAL REAGENTS  
**SAFETY DATA SHEET**

Date Prepared: September 2017  
Version No: 5

---

**1. IDENTIFICATION OF THE MATERIAL AND SUPPLIER**

---

Product Name: Ammonia Solution 10% v/v  
Product Code: 5483  
Other Names: Nil  
Uses: Analytical Reagent

Supplier: Australian Chemical Reagents  
38-50 Bedford Street Gillman SA 5013

Contacts: Telephone: 61 08 84402000  
Fax: 61 08 84402001  
Emergency Phone: 61 08 84402000 Mon – Fri 8:30am – 5:00pm

---

**2. HAZARDS INFORMATION**

---

**GHS Classification**

Skin Corrosion/Irritation: Category 2  
Serious Eye Damage/Irritation: Category 2A  
Specific Target Organ Toxicity – Single Exposure (Respiratory Tract) 3

**Signal Word(s)**  
**Pictogram(s)**

WARNING

**Hazard Statement(s)**

H315 Causes skin irritation.  
H319 Causes serious eye irritation.  
H335 May cause respiratory irritation

**Precautionary Statement(s)**  
**Preventative**

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.  
P264 Wash thoroughly after handling.  
P271 Use only outdoors or in a well-ventilated.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.

**Response**

P302+P352 IF ON SKIN: Wash with plenty of soap and water.  
P337+P313 If skin irritation persists: Get medical advice/attention.  
P362 Take off contaminated clothing and wash before reuse.  
P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P312 Call a POISON CENTER or doctor/physician if you feel unwell.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P337+P313 If eye irritation persists: Get medical advice/attention.

**Storage** P403+P233 Store in a well-ventilated place. Keep container tightly closed.  
P405 Store locked up.

**Disposal** P501 Dispose of contents/container to an approved waste disposal plant.

---

### 3. COMPOSITION / INFORMATION ON INGREDIENTS

---

#### Ingredients :

Chemical Entity	CAS No	Proportion
Ammonium hydroxide	[1336-21-6]	3%wt.
Water	[7732-18-5]	to 100%

---

### 4. FIRST AID MEASURES

---

Safety showers and eye wash facilities should be provided.

#### **Swallowed :**

If conscious wash out mouth with water. Seek medical advice. Show this SDS to medical practitioner.

#### **Eye :**

Immediately hold eyelids open and flood with water for at least 15 minutes. Obtain medical aid. Show this SDS to medical practitioner.

#### **Skin :**

Remove contaminated clothing. Immediately wash skin thoroughly with water and mild soap. Seek medical advice if irritation persists. Show this SDS to medical practitioner. Launder clothing before reuse.

#### **Inhaled :**

Remove from contaminated air. Maintain breathing with artificial respiration if necessary. Seek medical assistance. Show this SDS to a doctor.

---

### 5. FIRE FIGHTING MEASURES

---

#### **Suitable Extinguishing Media:**

Water spray. Carbon dioxide, dry chemical powder, or appropriate foam.

#### **Hazards From Combustion Products:**

Solutions will not burn or support combustion. Decomposition products include ammonia fumes.

#### **Precautions For Fire Fighters and Special Protective Equipment:**

Fire fighters and others who may be exposed to combustion products during fire should wear full protective clothing including positive pressure self-contained breathing apparatus (SCBA). Wear SCBA with full face-piece, operated in positive pressure mode when fighting fires.

---

### 6. ACCIDENTAL RELEASE MEASURES

---

#### **Emergency procedures:**

Prevent from entering waterways. Ventilate area. Restrict access to area. Remove chemicals that can react with the spilled material.

#### **Methods and materials for containment and clean up:**

Ventilate area. Use inert material such as sand or earth to contain spill or leak. Absorb spills with chemical absorber or vermiculite and dispose of in accordance with local regulations.

---

### 7. HANDLING AND STORAGE

---

#### **Precautions for Safe Handling:**

Do not get in eyes, on skin, on clothing. Avoid prolonged or repeated exposure.

#### **Conditions for Safe Storage:**

Store sealed in original container in a cool well ventilated situation away from foods and other chemicals. Do not store in direct sunlight. Observe good hygiene and housekeeping practices.

---

---

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

---

### National Exposure Standards:

Safe Work Australia – Ammonia TWA - 17 mg/m<sup>3</sup> STEL - 24 mg/m<sup>3</sup>

**Biological Limit Values:** No data available.

### Engineering Controls:

Do not use in confined spaces. Use only chemical fume hood.

### Personal Protective Equipment (PPE):

The use of nitrile or neoprene gloves complying with AS 2161 and the use of faceshield, chemical goggles or safety glasses with side shield protection complying with AS/NZS 1337 is recommended.

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

---

<b>Appearance :</b>	Clear liquid
<b>Odour:</b>	Ammonia
<b>pH:</b>	12
<b>Boiling Point (°C) :</b>	Not applicable
<b>Freezing/melting Point:</b>	Not applicable
<b>Vapour Pressure (mm of Hg @ 25°C) :</b>	Not applicable
<b>Vapour Density:</b>	Not applicable
<b>Specific Gravity :</b>	1
<b>Flash Point (°C) :</b>	Not flammable
<b>Flammability Limits (%) :</b>	Not flammable
<b>Solubility in Water (g/L) :</b>	Soluble

---

## 10. STABILITY AND REACTIVITY

---

### Chemical stability:

Stable.

### Conditions to avoid:

Heat.

### Incompatible materials:

Bases, iodine, strong acids, oxidizing agents

### Hazardous decomposition products:

Refer to section 5 (Fire Fighting Measures).

### Hazardous reactions:

Hazardous polymerization will not occur.

---

## 11. TOXICOLOGICAL INFORMATION

---

### Health Effects:

**Swallowed :** May cause burns. May be harmful if swallowed. Ammonium hydroxide orl-hmn LDLo 43mg/kg.

**Eye :** May be destructive to eye tissue. Ammonium hydroxide eye-rbt 750 ug sev. eye-rbt 100mg rinse sev.

**Skin :** May causes burn. May be harmful if absorbed through skin.

**Inhaled :** May be harmful if inhaled. May cause burning sensation, coughing, wheezing. Ammonium Hydroxide ihl-hmn LCLo 5000ppm

### Chronic Effects:

.

---

## 12. ECOLOGICAL INFORMATION

---

**Ecotoxicity:**

No data available. Avoid release to the environment.

**Persistence and degradability:**

No data available.

**Mobility:**

No data available.

---

## 13. DISPOSAL CONSIDERATIONS

---

Contact a licensed professional waste disposal service to dispose of this material. Observe all federal, state and local environmental regulations.

---

## 14. TRANSPORT INFORMATION

---

**UN Number:** None Allocated

**UN Proper Shipping Name:** None Allocated

**Class and subsidiary risk(s):** None Allocated

**Packing Group:** None Allocated

**Hazchem Code:** None Allocated

**Special precautions for user :** Nil

---

## 15. REGULATORY INFORMATION

---

**Standard for the Uniform Scheduling of Medicines and Poisons (SUSMP):**

Schedule 6

---

## 16. OTHER INFORMATION

---

**Disclaimer:**

All information given by the Company is offered in good faith and is believed to the best of our knowledge to be accurate. However this information is offered without warranty representation inducement or licence and the Company does not assume legal responsibility for reliance upon the same.

Every person dealing with the materials referred to herein does so at his or her own risk absolutely and must make independent determinations of suitability and completeness of information from all sources to ensure their proper use.