



**Product Specification**  
**Multi-purpose container, 120 ml, (LxØ): 105 x 44 mm, graduated, PP**



**Product description**

Order number	75.9922.420
Product description	Multi-purpose container, max. working volume: 120 ml, (LxØ): 105 x 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, sterile, 250 piece(s)/bag

**Product characteristics**

Graduation	yes
Closure type	screw cap
Base shape	flat base
Volume of work	120 ml

**Size**

Sample volume	120 ml
Diameter	44 mm
Height of product	105 mm
Length of product	105 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Multi-purpose container, 120 ml, (LxØ): 105 x 44 mm, graduated, PP**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	250
Piece(s) / pallet	4000
Case volume	0.0816 cbm
Weight of product	0.0188 kg
Weight of case	5.62 kg
EAN of inner box	4038917151863
EAN case	4038917038461

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2023-08-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**