



### Product description

Order number	72.991.002
Product description	PCR strip of 8, up to 200 µl, High Profile, PCR Performance Tested, transparent, material: PP, with anti-contamination shield, flat cap, cap attached, suitable for qPCR, 120 piece(s)/minigrip bag

### Product characteristics

Graduation	no
Cap	flat cap
Closure type	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	High Profile
Number of wells	8
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Sample volume	0.2 ml
Max. volume	200 µl
Diameter	7.5 mm
Width of product	10.2 mm
Height of product	20.65 mm
Length of product	72.4 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

### Packaging

Minimum order qty.	480
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	120
Piece (s) / outer case	480
Piece(s) / pallet	34560
Depth of case	358 mm
Width of case	323 mm
Height of case	140 mm
Case volume	0.0162 cbm
Weight of product	0.0015 kg
Weight of case	0.99 kg
EAN of inner box	4038917124294
EAN case	4038917067546

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

