



Product description

Order number	62.9924.283
Product description	Tube, working volume: 10 ml, (LxØ): 100 x 16 mm, material: PP, conical base, transparent, HD-PE, cap: natural, cap enclosed, with print, label/print: white, with graduations, 100 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Diameter	16 mm
Length of product	100 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Case volume	0.0551 cbm
Weight of product	0.0054 kg
Weight of case	6.2 kg
EAN of inner box	4038917155182
EAN case	4038917027915

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

